

corlife oHG

Feodor-Lynen-Str. 23
D-30625 Hannover

Phone: +49 (0)511 563 539 40

Fax: +49 (0)511 563 539 55

info@corlife.eu

www.corlife.eu

Commerzbank AG, Sarstedt

IBAN: DE91 2504 0066 0256 0001 00

BIC: COBADEFF138

USt-IdNr./VAT: DE247760300

Amtsgericht Hannover, HRA 202628

corlife oHG · Feodor-Lynen-Str. 23 · 30625 Hannover · Germany

Practitioners of patients having received a
“Decellularized Human Pulmonary Valve, Espoir PV” and/or a
“Decellularized Human Aortic Valve”

Compatibility with Diagnostic Imaging

October 17th, 2016


Dear doctor,

The “Decellularized Human Pulmonary Valve, Espoir PV” and the “Decellularized Human Aortic Valve, Arise AV” are allografts, which have been treated with detergents and buffers to remove cells of the donor. The implants consist predominantly of collagens and elastic fibres, the extracellular matrix.

There is no contraindication to diagnostic imaging, such as MRI, CT and ultrasound. However, the patient may have received other implants which are not compatible with imaging techniques.

Please do not hesitate to contact us if you have any questions.

Sincerely,



Dr. Michael Harder
geschäftsführender Gesellschafter